Chemistry Lab: Titration

**Introduction:** 

**Neutralization Reaction:** a reaction where an **acid** and a **base** react to form **water** and a **salt**; can be used to create a pure sample of a salt or to determine the concentration of an acid or a base in a solution.

Examples: HCl + NaOH  $\rightarrow$  NaCl + H<sub>2</sub>O H<sub>2</sub>SO<sub>4</sub> + 2 NaOH  $\rightarrow$  2 H<sub>2</sub>O + Na<sub>2</sub>SO<sub>4</sub>

**Equivalence point** (The point of neutralization in a titration): The amount of acid (or base) that will give one mole of hydrogen (or hydroxide) ions;

 $HCl + NaOH \rightarrow NaCl + H_2O$ 

In this example: one mole of HCl is one equivalent of HCl

Because moles of H+ from the acid are equal to moles of OH- from base, in the balanced equation, we can use the dilution equation to solve for unknown molarities.  $M_1V_1=M_2V_2$ 

<b>Balnced equation:</b>	HCl	+	NaOH	<b>→</b>	NaCl	+	H <sub>2</sub> O
Dilution equation:	M <sub>HCl</sub>	x V <sub>HCl</sub> =	M <sub>NaOH</sub> x V	NaOH	at		equivalnce point
Calculation:	? 2	x 35.3mL	.102M x 32	.8mL =	.102 x 32.	8 / 35	.5 = .0948 (3 sigfigs!)

**Objectives:** 1. Determine the M<sub>HCl</sub> for each titration

- 2. Average three molarities
- 3. Have me sign

Data:

Titration	HCl	NaOH	Molarity

1	Vol <sub>F</sub> HCl = Vol <sub>I</sub> HCl = Vol HCl used =	Vol NaOH <sub>F</sub> = Vol NaOH <sub>I</sub> = Vol NaOH used=	$M_1$ HCl =		
2	Vol <sub>F</sub> HCl = Vol <sub>I</sub> HCl = Vol HCl used =	Vol NaOH <sub>F</sub> = Vol NaOH <sub>I</sub> = Vol NaOH used=	M <sub>2</sub> HCl =		
3	Vol <sub>F</sub> HCl = Vol <sub>I</sub> HCl = Vol HCl used =	Vol NaOH <sub>F</sub> = Vol NaOH <sub>I</sub> = Vol NaOH used=	M <sub>3</sub> HCl =		
Average Molarity of HCl =					